



*Innovative Solutions in
Cryogenic Instrumentation*

USB 1.1 / 2.0 Option 4001-001



Features:

- Compatible with USB 1.1 / 2.0 equipment.
- Connect to the existing serial port on a Cryo-con instrument to add fully functional USB serial port emulation.
- Drivers provided for Windows, Macintosh and Linux computers. Simple setup and installation.
- High speed. Up to 57.6K baud.
- Uses existing LabView™ serial port drivers.
- No configuration required.
- Field installable at any time.

Contact Information

Cryogenic Control Systems, Inc.
PO Box 7012
Rancho Santa Fe, CA 92067
Tel: (858) 756-3900 Fax: (858) 759-3515
E-mail: sales@cryocon.com Web: www.cryocon.com

CalGen® and Cryo-Con® are registered trademarks of Cryogenic Control Systems, Inc. All other product and company names are trademarks or trade names of their respective companies.

©Cryogenic Control Systems, Inc. USBds 12/09
Specifications subject to change without notice